A J COPE & SON LTD

The Oval

London E2 9DU

Telephone: +44 (0) 20 7729 2405 Facsimile: +44 (0) 20 7729 2657

Internet: www.ajcope.co.uk

ELECTRONIC CATALOGUE CUSTOMER'S AMENDMENTS ISSUE 3 — March 2011

Page 1 of 7 pages

	G. IG
Amendment	
CANCELLATION	
PJ865-44 Cancelled	d. offer P1865-55
	, one i good ee
NEW ITEM	
Bench Meter H1221 ▲ pH. temperatu	13 re (°C) and mV ranges
▲ Simultaneous d	display of pH/mV and temperature
	or 3 point calibration with buffer recognition of up to 5 pH values
▲ Calibration ren	
▲ Manual or auto	omatic temperature compensation
Model	HI2213
Range	pH -2 to 16 (0.01)
(resolution)	mV -699.9 to 2000 (0.1)
	°C -20.0 to 120.0 (0.1)
Accuracy	pH ±0.01
	mV ±0.2 (-699.9 to 699.9 mV) ±1 (rest of range)
	$^{\circ}$ C ± 0.2
Calibration	Automatic 1, 2 or 3 point, pH4.01, 6.86, 7.01,
Canoration	9.18, 10.01 buffers recognised
Temp. compensation	
Input impedance	Ω 10 ¹²
Power	12V d.c. mains adapter
Overall L x W x D	mm 235 x 222 x 109
Bench Meter, HI-22	
	ous PJ865-44, but with wider and improved pH, mV and temperature ranges d now with 1, 2 or 3 point calibration.
PJ865-55 HI-	2213
Spares and accessor	ories
	ss bodied electrode
PJ862-22 Ten	mperature probe
CANCELLATION	
BB400 (2 C "	J. effert DD240-52
PP280-63 Cancelled	a, oner PP240-53

ISSUE 3 — March 2011 Page 2 of 7 pages

CANCELLATION

PR060-80 Cancelled, offer PR060-85

NEW ITEMS

Repeating Pipettor, HandyStep® S

Generally as previous PR060-80, but with a dispensed volume range from $2\mu l$ to $50000\mu l$ and nine separate volume settings.

- **▲** Ergonomic design gives comfortable handling.
- **▲** Nine separate volume settings
- ▲ Piston lock ensures only complete step volumes are dispensed
- ▲ Greatly reduced risk of pipettor contamination, or sample carry-over
- **▲** Factory calibrated
- **▲** Maintenance free

HandyStep® S Control Unit

As described. Supplied with wall support and instructions. Requires, but is not supplied with, PD-Tips for operation.

PR060-8	5	HandyStep® S Repeating Pipettor									
Control Unit setting	Max. no. of steps	Dispensing volume per step, µl									
	•	0.1ml	0.5ml	1ml	1.25ml	2.5ml	5ml	10ml	12.5ml	25ml	50ml
		PD-tip	PD-tip	PD-tip	PD-tip	PD-tip	PD-tip	PD-tip	PD-tip	PD-tip	PD-tip
		•	•	•	•	•	•	•	•	•	•
1	49	2	10	20	25	50	100	200	200	500	1000
1.5	32	3	15	30	37.5	75	150	300	375	750	1500
2	24	4	20	40	50	100	200	400	500	1000	2000
2.5	19	5	25	50	62.5	125	250	500	625	1250	2500
3	15	6	30	60	75	150	300	600	750	1500	3000
3.5	13	7	35	70	87.5	175	350	700	875	1750	3500
4	11	8	40	80	100	200	400	800	1000	2000	4000
4.5	10	9	45	90	112.5	225	450	900	1125	2250	4500
5	9	10	50	100	125	250	500	1000	1250	2500	5000
Accuracy	(A) %	±1.6	± 0.8	± 0.8	±0.8	± 0.7	±0.5	±0.4	±0.3	±0.3	±0.3
Variance	(V) %	≤5.0 to	≤1.4 to	≤1.0 to	≤0.8 to	≤0.8 to	≤0.6 to	≤0.3 to	≤0.4 to	≤0.4 to	≤0.4 to

PD-Tips, non sterile

≤2.0

≤0.6

≤0.4

As described. Supplied loose in packs as indicated. The 25ml and 50ml sizes require, and are supplied with, an autoclavable adapter.

≤0.2

≤0.2

≤0.2

≤0.2

≤0.2

≤0.15

≤0.3

	Size	Pack
	ml	qty.
PR064-08	0.1	100
PR064-10	0.5	100
PR064-12	1	100
PR064-14	1.25	100
PR064-18	2.5	100
PR064-22	5	100
PR064-24	10	100
PR064-26	12.5	100
PR064-30	25	50
PR064-34	50	25

ISSUE 3 — March 2011 Page 3 of 7 pages

PD-Tips, sterile

As described. Endotoxin free and beta irradiated in accordance with ISO 8655. The sterility meets USP XXIII requirements. Individually wrapped within outer packs as indicated. The 25ml and 50ml sizes require, and are supplied with, an autoclavable adaptor.

	Size ml	Pack qty.
PR067-08	0.1	100
PR067-10	0.5	100
PR067-12	1	100
PR067-14	1.25	100
PR067-18	2.5	100
PR067-22	5	100
PR067-24	10	100
PR067-26	12.5	100
PR067-30	25	50
PR067-34	50	25

AMENDMENT

SA300-60 Cancelled, offer SA300-62

NEW ITEM

SA300-62 Hoods, TYVEK, white, one-size, pack of 20

CANCELLATIONS

SB896-10 SB896-20 SB896-30	Cancelled, offer SB896-15 Cancelled, offer SB896-25
SB890-30	Cancelled, offer SB896-35
SB897-30	Cancelled, offer SB897-35
SB897-40	Cancelled, offer SB897-45
SB897-50	Cancelled, offer SB897-55
SB898-40	Cancelled, offer SB898-45
SB898-50	Cancelled, offer SB898-55

NEW ITEMS

Stainless Steel Safety Cans

For carrying, dispensing and storage of corrosive, flammable or contamination-sensitive liquids.

- ▲ Tüv/GS* approved
- ▲ Grade 1.4571 stainless steel construction with powder-coated, ergonomic, aluminium handle
- ▲ Leakproof closure with Teflon gasket
- ▲ Spring-operated pressure/vacuum relief valve
- ▲ Perforated stainless steel flame arrester in spout resists flashback ignition
- **▲** Base formed in one-piece
- ▲ Neoprene bumper guard
- **▲** Leakproof closure with Teflon gasket
- ▲ Fitted with carry handle

ISSUE 3 — March 2011 Page 4 of 7 pages

Safety Cans

Stainless steel construction with choice of self-closing, leakproof dispenser cap and nozzle which can be removed to allow filling or $32 \text{mm} (1^{1}/4^{"})$ screw cap.

	Capacity	Overall Ref
	litres	dia. x H, mm
Self-closing,	dispenser cap	
SB896-15	i	111 x 235
SB896-25	2	131 x 295
SB896-35	5	165 x 400
Screw cap		
SB896-40	1	111 x 235
SB896-50	2	131 x 295
SB896-60	5	165 x 400

Safety Canisters

Stainless steel construction with leakproof 37mm $(1^{1}/_{2}")$ screw closure, stainless steel flame arrester, carrying handle and neoprene bumper guard.

	Capacity	Overall Ref
	litres	W x D x H, mm
SB897-35	5	130 x 280 x 310
SB897-45	10	130 x 280 x 380
SB897-55	20	175 x 345 x 470

Dispensing Canisters

Generally as SB897 series but with self-closing, leakproof tap and separate vent to allow dispensing from a shelf or an accessory dispensing station.

	Capacity	Overall		
	litres	W x D x H, mm		
SB898-45	10	130 x 280 x 380		
SB898-55	20	175 x 345 x 470		

Stainless Steel Filling Funnels

For use with SB896 to SB898 canisters and cans, unless indicated.

SB898-60	Straight, plug-fit
SB898-62	Angled, plug-fit
	0 /1 0

SB898-64 Angled, with screw connection and overfill protection for SB895/SB898-series only

Accessory Dispensing Station

For use with up to two SB898 series dispensing cans. Free standing or wall mounting steel frame with removable stainless steel drip tray. Overall dimensions $555 \times 405 \times 600$ mm (W x D x H). Without cans.

SB898-95 Dispensing station

CANCELLATION

SC360-58 Cancelled, offer SC360-53

CANCELLATION

SD057-29 Cancelled, no replacement

ISSUE 3 — March 2011 Page 5 of 7 pages

AMENDMENTS

SM805-10 to -18 series is now grey, not white as stated in the catalogue.

All other details are unchanged.

CANCELLATIONS

SV340-15 Cancelled, offer SV340-20

SV358-10 Cancelled, no replacement SV358-15 Cancelled, no replacement

REINSTATEMENT

Spare PTFE keys for Pyrex stopcocks

Bore Type

mm

SV358-20 4 T-bore

CANCELLATIONS

TJ847-50 Cancelled, offer TJ847-60 TJ847-55 Cancelled, offer TJ847-65 TJ847-70 Cancelled, offer TJ847-80

NEW ITEMS

Dataloggers, 175, 176

Long term, reliable temperature loggers with robust ABS casing providing either IP65 or IP68 protection, push button start, SD card or USB data connection and flexible data handling options.

- Model 175 T1 Single channel, with internal NTC sensor and LCD readout
- Model 175 T2 Two channel, with internal NTC sensor, external NTC probe socket for use with accessory probes and twin-line LCD readout
- Model 176 T1 Single channel, expanded memory version, with internal sensor but without display
- ▲ Internal interfaces for both USB or SD card data collection
- ▲ Measuring time cycle adjustable from 10 seconds to 24 hours (175 series) or 1 second to 24 hours (model 176)

Model			175 T1	175 T2	176 T1
Ranges	internal sensor	$^{\circ}\mathbf{C}$	-35 to +55	-35 to +55	-35 to $+70$
	external sensor	$^{\circ}\mathbf{C}$	_	40 to +120	_
Resolution		°C	0.1 — both mo	odels	0.01
Memory	readings		1 million - bot	h models	2 million
Memory Overall, W x 1	0		1 million - bot 53 x 27 x 89	h models 53 x 27 x 89	2 million 63 x 33 x 103
•	0	g			

ISSUE 3 — March 2011 Page 6 of 7 pages

As described. Supplied with both USB and SD card data collection interfaces, wall holder, batteries and instructions. Requires, but is not supplied with, either an accessory USB cable or SD card for data transfer.

TJ847-60 Testo 175 T1 TJ847-65 Testo 175 T2 TJ847-80 Testo 176 T1

Accessories

Accessory probes for TJ847-65

With 1.5m lead. Dimensions are probe length x diameter.

TJ802-14 Immersion, 110 x 4mm TJ802-17 Air, 110 x 4mm

TJ812-20 Laboratory, glass sheathed, 200 x 6mm

TJ847-87 USB cable for connecting 175 and 176

dataloggers to the user's PC.

Comsoft 4 professional software

With data handling, database, graphical and trend analysis functions.

TJ847-98 Comsoft 4

T H E E N D

Page 7 of 7 Pages

March 2011

CAT. NO.	PRICE £	Per	C	ODE
PJ865-55	458.00	each		FB
PR060-85	190.00	each		OG
PR064-08	66.50	pack		AO
PR064-12	60.00	pack		AO
PR064-24	60.00	pack		AO
PR067-08	132.00	pack		AO
PR067-12	118.00	pack		AO
PR067-24	118.00	pack		AO
SA300-62	39.80	pack		ЕО
SB896-15	167.00	each		ЕО
SB896-25	180.00	each		ЕО
SB896-35	194.00	each		EO
SB896-40	129.00	each		ЕО
SB896-50	148.00	each		EO
SB896-60	170.00	each		EO
SB897-35	199.00	each		EO
SB897-45	235.00	each		EO
SB897-55	280.00	each		EO
SB898-45	290.00	each		EO
SB898-55	320.00	each		EO
SB898-60	52.00	each		EO
SB898-62	58.00	each		EO
SB898-64	103.00	each		EO
SV358-20	14.00	each	7.00 each per 5	ОН
TJ847-60	110.00	each		RO
TJ847-65	135.00	each		RO
TJ847-80	299.00	each		RO
TJ847-87	20.90	each		RO
TJ847-98	272.00	each		RO

T H E E N D