

ELECTRONIC CATALOGUE
CUSTOMER'S AMENDMENTS
ISSUE 9 — September 2012
Page 1 of 18 pages

| | <i>Amendment</i> |
|--|---|
| | <p>CANCELLATIONS</p> <p>AD180-06 Cancelled, no replacement AD190-10 Cancelled, no replacement AD195-18 Cancelled, no replacement</p> <hr/> |
| | <p>CANCELLATION</p> <p>BN240-16 Cancelled, no replacement</p> <hr/> |
| | <p>CANCELLATION</p> <p>BR154-42 Cancelled, offer BR154-45</p> <hr/> |
| | <p>CANCELLATION</p> <p>BR406-48 Cancelled, offer BR406-50</p> <p>REINSTATEMENT</p> <p><u>Bottles, Wide mouth Powder, Amber glass</u> With white polypropylene wadless cap.</p> <p>BR406-50 Bottles, Wide mouth powder, amber glass, 500ml, R3/58 neck, pack of 28</p> |

CANCELLATION

BR756-60 Cancelled, offer BR756-62

NEW ITEM

Bottle, polyethylene, leakproof

High quality, sturdy translucent bottle, designed for laboratory use at an economical price. Wide mouth with square shoulders and leakproof cap. Supplied in pack as indicated.

| | Cap. ml | Pack qty. |
|-----------------|--------------------|----------------------|
| BR756-62 | 2000 | 1 |

CANCELLATION

BS800-10 Cancelled, no replacement

NEW ITEMS

Safety Caps

SCAT Europe

For safe closure and venting of supply reservoirs, particularly when handling solvents for HPLC.

- Replaceable air inlet filter stops solvent contamination from dust**
- Halts solvent evaporation, reducing cost and contamination of the laboratory atmosphere by fumes**
- Precision adapters allow insertion without crimping and fixed positioning of capillary tubing**
- Easy reservoir exchange**
- Manufactured in PTFE and HDPE for resistance to aggressive solvents**
- Caps can be unscrewed without tubing removal**
- Fire-resistant versions available**

As described. Supplied singly, with one inlet filter (6-monthly replacement recommended) and capillary tubing connectors for tubing with outside diameters as indicated.

| | Cap type | Connectors for tubing |
|------------------------|-------------|-------------------------------|
| <u>Standard</u> | | |
| CK780-04 | GL45 | 1 x 3.2mm o.d. |
| CK780-08 | GL45 | 1 x 3.2mm o.d. with stopcock |
| CK780-12 | GL45 | 2 x 3.2mm o.d. |
| CK780-16 | GL45 | 3 x 3.2mm o.d. |
| CK780-20 | GL45 | 4 x 3.2mm o.d. |
| CK780-24 | GL45 | 4 x 3.2mm o.d. with stopcocks |
| CK780-28 | GL45 | 6 x 3.2mm o.d. |

| | Cap type | Connectors for tubing |
|------------------------------|-------------|-------------------------------|
| <u>Fire-resistant</u> | | |
| CK784-04 | GL45 | 1 x 3.2mm o.d. |
| CK784-08 | GL45 | 1 x 3.2mm o.d. with stopcock |
| CK784-14 | GL45 | 2 x 3.2mm o.d. with stopcocks |
| CK784-24 | GL45 | 4 x 3.2mm o.d. with stopcocks |
| CK784-28 | GL45 | 6 x 3.2mm o.d. with stopcocks |

Accessory inlet filters

| | |
|----------|---------------------------|
| CK785-10 | Standard, pack of 1 |
| CK785-20 | Standard, pack of 10 |
| CK785-40 | Fire-resistant, pack of 1 |

Safety Waste Caps

SCAT Europe

For safe closure and venting of liquid waste reservoirs.

- Manufactured in PTFE and HDPE for resistance to aggressive solvents, acids and bases
- Easy reservoir exchange
- Caps can be unscrewed without tubing removal
- Precision adapters allow insertion without crimping and fixed positioning of capillary tubing
- Accessory, activated charcoal, exhaust filters absorb up to 99% of volatiles and are available in 3 sizes to accommodate alternative sizes of reservoir

As described. Supplied singly with connectors for one exhaust filter, capillary tubing with outside diameters and tubing ferrule as indicated. Without exhaust filter.

| | Cap type | Connectors for tubing |
|----------|----------|---|
| CK790-03 | GL45 | 3 x 2.3/3.2mm o.d. |
| CK790-06 | GL45 | 2 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |
| CK790-07 | GL45* | 2 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |
| CK790-09 | S51 | 2 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |
| CK790-12 | 53B | 4 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm o.d. |
| CK790-15 | S55 | 2 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |
| CK790-18 | S60/61 | 3 x 2.3/3.2mm o.d. |
| CK790-19 | S60/61 | 2 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |
| CK790-22 | 83B | 4 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |
| CK790-26 | S90 | 4 x 2.3/3.2mm o.d. 1 x L-ferrule 6.4mm i.d. |

*Fire-resistant

Accessory Exhaust filters

With activated charcoal absorbent.

| | Size | Filling, g | For vessels | Pack |
|----------|------|------------|-------------|------|
| CK794-05 | S | 24 | <5L | 1 |
| CK794-07 | S | 24 | <5L | 4 |
| CK794-12 | M | 48 | <20L | 1 |
| CK794-16 | M | 48 | >20L | 3 |
| CK794-19 | M* | 48 | >20L | 1 |
| CK794-24 | L | 100 | >20L | 1 |

*Fire-resistant

Multicollector Manifolds

Allow multiple waste streams to be routed through one inlet in the waste cap.

| | |
|----------|--|
| CK796-10 | 2 x 3.2mm o.d., 1 x 6.4mm i.d. side inlets |
| CK796-15 | 3 x 3.2mm o.d., side inlets |
| CK796-20 | 2 x 6.4mm i.d., side inlets |
| CK796-25 | 3 x 6.4mm i.d., side inlets |
| CK796-55 | 8 x 3.2mm o.d., top inlets |

Fittings and Ferrules for Safety and Safety Waste Caps and Manifolds

| | |
|----------|---|
| CK798-02 | Fittings, PP, for 1.6mm, 2.3mm and 3.2mm o.d. capillary tubes, pack of 10 |
| CK798-04 | Fittings with ferrule, PTFE, for 1.6mm o.d. capillary tubing, pack of 10 |
| CK798-06 | Fittings with ferrule, PTFE, for 2.3mm o.d. capillary tubing, pack of 10 |
| CK798-08 | Fittings with ferrule, PTFE, for 3.2mm o.d. capillary tubing, pack of 10 |

- CK798-10** Fitting with ferrule, PTFE, for 4.76mm (3/16") o.d. capillary tubing, pack of 1
CK798-22 Plugs, PTFE, for 3.2mm i.d. fittings, pack of 50
CK798-24 Plugs, PTFE, for 4.76mm (3/16") fittings, pack of 10
CK798-35 Ferrules, PTFE, for 1.6mm i.d. fittings, pack of 10
CK798-37 Ferrules, PTFE, for 2.3mm i.d. fittings, pack of 10
CK798-39 Ferrules, PTFE, for 3.2mm i.d. fittings, pack of 10

Capillary Tubing

PTFE. In lengths of 3 x 3m.

- CK799-10** 1.6mm o.d., 1.0mm i.d.
CK799-20 2.3mm o.d., 1.7mm i.d.
CK799-30 3.2mm o.d., 1.6mm i.d.
CK799-65 Inlet filters for HPLC solvents, PFA/PTFE to fit 3.2mm o.d. tubing, pack of 5

CANCELLATIONS

- CX304-23** Cancelled, offer **CX304-26**
CX308-35 Cancelled, offer **CX308-32**

CANCELLATIONS

- DA450-30** Cancelled, offer **DA450-40**
DA454-15 Cancelled, now included with **DA450-40**

NEW ITEMS

Gerber Centrifuge, Universal

Generally as previous **DA450-30** but with speed range 200 to 1300rpm and supplied with 12-place universal swing-out rotor. For 230V 50/60Hz single phase supplies. Without buckets or heating module.

- DA450-40** Universal Centrifuge

CANCELLATION

DT330-10 Cancelled, offer DT330-15

NEW ITEM

Dissolved Oxygen Meter, Model HI-9146

Generally as previous DT330-10, but with backlit display feature, altitude and salinity compensation, on-screen user guidance and GLP features.

| | | |
|---------------------------|--------------------------|-----------------------------------|
| Catalogue No. | | DT330-15 |
| Model | | HI-9146 |
| Range x resolution | ppmO₂ | 0 to 45 x 0.01 |
| | %O₂ | 0 to 300 x 0.1 |
| | temperature °C | 0 to +50 x 0.1 |
| Accuracy | ppm%O₂ | ±1.5% of full scale |
| | temperature °C | ±0.5 |
| Calibration | | Automatic in saturated air |
| Power | | 3 x 1.5V AAA batteries |
| Dimensions | mm | 185 x 72 x 36 |
| Weight | g | 300 |

DO₂ Meter, Hanna

As described. With rugged carrying case, dissolved oxygen probe with 4m cable, batteries, spare membrane kit and 30ml electrolyte solution.

DT330-15 HI-9146

Spares

- DT335-15 DO₂ probe with 4m cable**
- DT335-20 Replacement membranes, pack of 5**
- DT335-24 Electrolyte solution, pack of 30ml**
- BL610-08 1.5V AAA battery (3 required)**

CANCELLATION

DT980-20 Cancelled, no replacement

CANCELLATION

FK830-35 Cancelled, no replacement

CANCELLATIONS

FT250-12 Cancelled, offer FT252-12
FT250-50 Cancelled, offer FT280-55

NEW ITEMS

Funnels, Conical, Polypropylene

General purpose, conical funnels. Autoclavable at 121°C and 1.1bar for 30 minutes. Dimensions given are top diameter (including flange) x maximum stem length x diameter. Supplied singly.

FT252-12 40 x 45 x 6mm
FT252-25 58 x 48 x 8mm
FT252-30 70 x 62 x 8mm
FT252-35 82 x 75 x 10mm
FT252-40 110 x 95 x 10mm
FT252-45 150 x 110 x 18mm

NEW ITEM

Silica Gel

Self indicating. Supplied in a pack of 25kg.

GG150-97 Silica gel

NEW ITEM

Gas Meters

Spare battery for digital counter on GG500-10 and GG500-15.

GG500-95 Spare battery

CANCELLATION

ME884-10 Cancelled, no replacement

CANCELLATION

ML700-15 Cancelled, offer ML700-10

AMENDMENTS

PD050-40 and -45 are now supplied in packs of 700.

All other details are unchanged.

CANCELLATIONS

PY250-10 Cancelled, offer PY260-10

PY250-30 Cancelled, offer PY260-30

PY250-45 Cancelled, no replacement

PY252-06 Cancelled, no replacement

PY252-08 Cancelled, no replacement

PY252-15 Cancelled, no replacement

NEW ITEMS

Scroll Pumps

Edwards

Innovative, portable, dry pumps designed to provide vacuum totally free of lubricants. Ultimate vacuum 0.007, 0.02 or 0.03mbar, depending on model selected

- No oils or greases of any kind required, ideal for “clean” processes**
- Meshing scroll mechanism**
- Patented bearing shield isolates bearings from exposure to the vacuum environment**
- Low noise and vibration levels**
- Very low maintenance**

| Catalogue No. | | PY260-10 | PY260-30 | PY260-45 | PY260-60 |
|-----------------------------------|-------------------|------------------------------|-----------|-----------|-----------|
| Model No. | | nXDS6i | nXDS10i | nXDS15i | nXDS20i |
| Ultimate vacuum | mbar | 0.02 | 0.007 | 0.007 | 0.03 |
| (total pressure) torr | | 0.015 | 0.005 | 0.005 | 0.022 |
| Displacement | m ³ /h | 6.8 | 12.7 | 17.1 | 28.0 |
| Maximum water vapour pumping rate | kg/h | 0.11 | 0.14 | 0.28 | 0.22 |
| Inlet/outlet connection | | NW25/NW25 | NW25/NW25 | NW25/NW25 | NW25/NW25 |
| Gas ballast on/off | | Yes | Yes | Yes | Yes |
| Noise level (at 1 metre) | db(A) | 52 | 52 | 52 | 52 |
| Vibration | mm/sec.(rms) | <4.5 | <4.5 | <4.5 | <4.5 |
| Overall, H x W x D | mm | 302 x 265 x 432 — all models | | | |
| Weight | kg | 26.2 | 25.8 | 25.2 | 26.6 |

Scroll Pumps, nXDS series

As described. With direct drive motor. For 100-127V and 200-240V 50/60Hz single phase supplies. Supplied without mains lead.

PY260-10 nXDS6i
PY260-30 nXDS10i
PY260-45 nXDS15i
PY260-60 nXDS20i

Accessories and Spares

PY262-06 Tip seal service kit for PY260-series

PY262-15 Silencer for PY260-series

Accessory mains leads

With IEX60320 connector and plug as indicated for use with Edwards pumps, PY260/PY330-series.

PY334-12 Without plug

PY334-14 With UK plug

PY334-18 With European Schuko plug

Clamping Ring

Stainless steel. For connecting gauges, adapters and pipe flanges. Requires PY787-06 centring ring seal.

PY785-06 NW25

Centring Ring Seal

Polymer, with nitrile O-ring. For use with PY785-06 clamping ring.

PY787-06 NW25

Adapters

Aluminium. NW25 connection.

PY789-10 Radiused elbow, 90°

PY789-15 T-Piece, 3-way

PY789-20 Cross, 4-way, 90° between limbs

PY789-40 Nozzle, for 12mm i.d. tubing

CANCELLATION

QSTQ Cancelled, no replacement

AMENDMENT

RE115-20 is now range 0 to 54% sugar (Brix) not 0 to 56% as stated in the catalogue.

All other details remain unchanged.

CANCELLATION

RE490-20 Cancelled, offer RE500-15

CANCELLATIONS

RG850-03 Cancelled, no replacement
RG850-06 Cancelled, no replacement
RG850-09 Cancelled, no replacement
RG850-12 Cancelled, no replacement
RG850-15 Cancelled, no replacement
RG850-18 Cancelled, no replacement
RG850-21 Cancelled, no replacement
RG850-24 Cancelled, no replacement
RG850-27 Cancelled, no replacement
RG850-30 Cancelled, no replacement
RG850-33 Cancelled, no replacement
RG850-36 Cancelled, no replacement
RG850-39 Cancelled, no replacement
RG850-42 Cancelled, no replacement
RG850-45 Cancelled, no replacement
RG850-48 Cancelled, no replacement
RG850-51 Cancelled, no replacement
RG850-54 Cancelled, no replacement
RG850-57 Cancelled, no replacement
RG850-61 Cancelled, no replacement
RG850-64 Cancelled, no replacement
RG850-67 Cancelled, no replacement
RG850-70 Cancelled, no replacement
RG850-73 Cancelled, no replacement
RG850-76 Cancelled, no replacement
RG850-79 Cancelled, no replacement
RG850-82 Cancelled, no replacement
RG850-85 Cancelled, no replacement
RG850-88 Cancelled, no replacement
RG850-91 Cancelled, no replacement
RG850-99 Cancelled, no replacement

=====

CANCELLATIONS

SE162-16 Cancelled, no replacement
SE162-19 Cancelled, no replacement
SE162-30 Cancelled, no replacement
SE162-33 Cancelled, no replacement
SE162-36 Cancelled, no replacement
SE162-39 Cancelled, no replacement

CANCELLATIONS

SJ379-10 Cancelled, offer SJ379-15
SJ379-45 Cancelled, no replacement
SJ379-54 Cancelled, no replacement
SJ379-63 Cancelled, no replacement
SJ379-72 Cancelled, no replacement
SJ379-81 Cancelled, no replacement
SJ379-90 Cancelled, no replacement
SJ379-95 Cancelled, no replacement

NEW ITEMS

Spectrophotometer, Double Beam Model 6850

Jenway

A true, high accuracy, double beam UV/Visible range instrument with variable 0.5, 1, 2, 4 and 5nm spectral bandwidths, using either the integrated controller or remotely by the user's PC via Windows® compatible Prism PC software supplied.

- ❑ **Six measurement modes:**
 - **Photometrics and multi-wavelength spectrum scanning (up to 10 wavelength can be measured, or 20 when using Prism PC software, with up to 0.1nm resolution and 2000nm/min. scan speed)**
 - **Spectrum scanning (up to 0.1nm resolution and 2000nm/min. scan speed)**
 - **Kinetics (up to 12 hours with time intervals of 0.1, 0.2, 0.5, 1, 2, 5, 10 or 30 seconds)**
 - **Quantitation (up to 3 wavelength points using up to 10 standards)**
 - **DNA/RNA (ratio, concentration, A320 correction)**
 - **Protein analysis (autopeaks and valleys, zoom, addition, subtraction, peak ratios, smoothing, area under curve, wavelength table, derivatives and overlay with PC software)**
- ❑ **Range: 190 to 1100nm**
- ❑ **Prism PC software is designed to give a comprehensive range of measurement modes whilst remaining intuitive, user-friendly and providing extensive post-measurement tools and functions**
- ❑ **Large graphical display**
- ❑ **Conforms to European Pharmacopoeia requirements**
- ❑ **Range of interchangeable sample chambers allows rapid reconfiguration of upgrade**

| | | |
|--------------------------------------|------------------|--|
| Catalogue No. | | SJ379-15 |
| Model | | 6850 |
| Wavelength range x resolution | nm | 190 to 1100 x 0.1 |
| Bandwidths | nm | 0.5, 1, 2, 4, 5 |
| Wavelength accuracy | nm | ±0.3 (at 0.5 and 1nm), ±0.5 (at 2, 4 and 5nm) |
| Straylight | % | <0.05% (at 220nm and 360nm) |
| Source switching range | nm | 325 to 370 |
| Photometric range | A | -0.3 to 3.0 |
| Photometric accuracy | A | ±0.002 |
| Photometric reproducibility | A | ±0.001 |
| Baseline stability | A/h | ±0.001 |
| Noise level | A | 0.0005 |
| Scan speed | nm/minute | 3 selectable, from 100 to 2000 |
| Interface | | USB and parallel port |
| Light source | | Tungsten-halogen/Deuterium |
| Supply requirements | | 220V 50/60Hz single phase supply |
| Overall | mm | 600 x 450 x 200 |
| Weight | kg | 22 |

As described. Supplied with single 10 x 10mm cuvette holder (sample and reference positions), instruction manual, Prism PC software with USB cable, 2 quartz and 4 glass cuvettes and dust cover. For 220V 50/60Hz single phase supplies.

SJ379-15 Model 6850, double beam, UV/Visible

Accessory and spare sample chambers

- SJ379-40 10 x 10mm path length cuvette holder**
- SJ379-50 Water heated 10 x 10mm single cell holder**
- SJ379-55 10 to 100mm adjustable path length cuvette holder**
- SJ379-60 Micro-cuvette holder**
- SJ379-65 8-position automatic cell changer**
- SJ379-85 Spare Tungsten-halogen lamp**
- SJ379-87 Spare Deuterium lamp**

CANCELLATION

SR875-40 Cancelled, special order only

CANCELLATIONS

SR946-10 Cancelled, special order only
SR946-20 Cancelled, special order only
SR946-30 Cancelled, special order only
SR946-40 Cancelled, special order only
SR946-75 Cancelled, special order only

CANCELLATION

SS365-20 Cancelled, offer SS345-49

CANCELLATION

SS600-20 Cancelled, offer SS705-10

CANCELLATIONS

SX235-02 Cancelled, no replacement
SX235-04 Cancelled, no replacement
SX235-06 Cancelled, no replacement
SX235-08 Cancelled, no replacement
SX235-10 Cancelled, no replacement
SX235-12 Cancelled, no replacement
SX235-14 Cancelled, no replacement
SX235-16 Cancelled, no replacement
SX235-18 Cancelled, no replacement
SX235-95 Cancelled, no replacement
SX235-96 Cancelled, no replacement
SX237-03 Cancelled, no replacement
SX237-05 Cancelled, no replacement
SX237-07 Cancelled, no replacement
SX237-09 Cancelled, no replacement
SX237-11 Cancelled, no replacement
SX237-13 Cancelled, no replacement
SX237-15 Cancelled, no replacement
SX237-17 Cancelled, no replacement
SX237-19 Cancelled, no replacement

NEW ITEM

Test Paper Books, Johnsons

Generally as TD300-20 but CE-marked in compliance with the IVD Directive.
Supplied in a pack of 10 books.

TD300-21 pH 1 to 11 universal, CE-marked

CANCELLATION

TG350-35 Cancelled, offer TG360-35

AMENDMENTS

TG300-35 and TG304-35 are now total immersion thermometers not partial immersion as stated in the catalogue.

All other details remain unchanged.

NEW ITEM

Multiwell Plates, Corning

U-bottom well without lid. Untreated polypropylene

Solvent resistant, non-sterile, untreated polypropylene, DMSO-, ethanol- and methanol-resistant. 96-well. Without lid. Sleeve quantity 25, pack quantity 100.

TP680-40 3365

CANCELLATIONS

TT817-04 Cancelled, no replacement
TT817-07 Cancelled, no replacement
TT817-11 Cancelled, no replacement
TT817-14 Cancelled, no replacement
TT817-75 Cancelled, no replacement
TT817-85 Cancelled, no replacement
TT817-95 Cancelled, no replacement
TT820-05 Cancelled, no replacement
TT820-15 Cancelled, no replacement
TT820-25 Cancelled, no replacement
TT820-50 Cancelled, no replacement
TT825-04 Cancelled, no replacement
TT825-09 Cancelled, no replacement

TT825-13 Cancelled, no replacement
TT825-17 Cancelled, no replacement
TT825-21 Cancelled, no replacement
TT825-25 Cancelled, no replacement
TT825-60 Cancelled, no replacement
TT830-10 Cancelled, no replacement
TT830-14 Cancelled, no replacement
TT830-18 Cancelled, no replacement
TT830-22 Cancelled, no replacement
TT830-26 Cancelled, no replacement
TT830-30 Cancelled, no replacement
TT830-55 Cancelled, no replacement

CANCELLATION

WD215-15 Cancelled, no replacement

NEW ITEM

WD215-27 Disposable bin liners for WD215-25, pack of 100

T H E E N D

PRICES FOR NEW ITEMS

Page 16 of 18 Pages

September 2012

| CAT. NO. | PRICE £ | Per | CODE |
|----------|---------|------|------|
| BR406-50 | 24.50 | pack | RO |
| BR756-62 | 3.35 | each | AE |
| CK780-04 | 89.00 | each | OO |
| CK780-08 | 131.60 | each | OO |
| CK780-12 | 100.50 | each | OO |
| CK780-16 | 106.40 | each | OO |
| CK780-20 | 142.40 | each | OO |
| CK780-24 | 201.50 | each | OO |
| CK780-28 | 143.75 | each | OO |
| CK784-04 | 145.60 | each | OO |
| CK784-08 | 205.00 | each | OO |
| CK784-14 | 249.00 | each | OO |
| CK784-24 | 269.00 | each | OO |
| CK784-28 | 357.00 | each | OO |
| CK785-10 | 27.00 | each | OO |
| CK785-20 | 264.00 | pack | OO |
| CK785-40 | 44.50 | each | OO |
| CK790-03 | 90.90 | each | OO |
| CK790-06 | 84.90 | each | OO |
| CK790-07 | 137.50 | each | OO |
| CK790-09 | 143.50 | each | OO |
| CK790-12 | 233.00 | each | OO |
| CK790-15 | 131.50 | each | OO |
| CK790-18 | 125.50 | each | OO |
| CK790-19 | 131.50 | each | OO |
| CK790-22 | 307.00 | each | OO |
| CK790-26 | 207.00 | each | OO |
| CK794-05 | 78.10 | each | OO |
| CK794-07 | 271.00 | pack | OO |
| CK794-12 | 112.00 | each | OO |
| CK794-16 | 292.00 | pack | OO |
| CK794-19 | 215.00 | each | OO |
| CK794-24 | 152.00 | each | OO |
| CK796-10 | 71.90 | each | OO |
| CK796-15 | 71.90 | each | OO |
| CK796-20 | 71.90 | each | OO |
| CK796-25 | 71.90 | each | OO |
| CK796-55 | 139.50 | each | OO |
| CK798-02 | 35.30 | pack | OO |
| CK798-04 | 84.90 | pack | OO |

PRICES FOR NEW ITEMS

Page 17 of 18 pages

| | | | |
|-----------------|----------------|-------------------|-----------|
| CK798-06 | 85.00 | pack | OO |
| CK798-08 | 85.00 | pack | OO |
| CK798-10 | 39.30 | each | OO |
| CK798-22 | 209.00 | pack | OO |
| CK798-24 | 87.70 | pack | OO |
| CK798-35 | 37.10 | pack | OO |
| CK798-37 | 37.10 | pack | OO |
| CK798-39 | 37.10 | pack | OO |
| CK799-10 | 11.65 | each | OO |
| CK799-20 | 11.65 | each | OO |
| CK799-30 | 16.45 | each | OO |
| CK799-65 | 62.40 | pack | OO |
| DA450-40 | 3165.00 | each | BO |
| DT330-15 | 457.00 | each | FB |
| FT252-12 | 1.00 | each | AO |
| | 0.35 | each (10+) | |
| FT252-25 | 1.00 | each | AO |
| | 0.40 | each (10+) | |
| FT252-30 | 1.00 | each | AO |
| | 0.75 | each (10+) | |
| FT252-35 | 1.00 | each | AO |
| | 0.75 | each (10+) | |
| FT252-40 | 1.80 | each | AO |
| | 1.50 | each (10+) | |
| FT252-45 | 4.55 | each | AO |
| | 3.80 | each (10+) | |
| GG150-97 | 178.50 | pack | OO |
| GG500-95 | 14.20 | each | BO |
| PY260-10 | 3195.00 | each | QA |
| PY260-30 | 3995.00 | each | QA |
| PY260-45 | 4795.00 | each | QA |
| PY260-60 | 5395.00 | each | QA |
| PY262-06 | 199.00 | each | QA |
| PY262-15 | 186.00 | each | QA |
| SJ379-15 | 6995.00 | each | KJ |
| SJ379-40 | 250.00 | each | KJ |
| SJ379-50 | 500.00 | each | KJ |
| SJ379-55 | 250.00 | each | KJ |
| SJ379-60 | 750.00 | each | KJ |
| SJ379-65 | 1000.00 | each | KJ |
| SJ379-85 | 100.00 | each | KJ |

PRICES FOR NEW ITEMS
Page 18 of 18 pages

| | | | |
|---------------------------|---------------|-------------------|-----------|
| SJ379-87 | 450.00 | each | KJ |
| TD300-21 | 6.05 | pack | RO |
| | 5.50 | pack (10+) | |
| TP680-40 | 173.30 | pack | OP |
| WD215-27 | 7.00 | pack | RO |
| | 6.60 | pack (10+) | |
| <i>T H E E N D</i> | | | |