A J COPE & SON LTD

The Oval

**London E2 9DU** 

Telephone: +44 (0) 20 7729 2405 Facsimile: +44 (0) 20 7729 2657

**Internet:** www.ajcope.co.uk

# ELECTRONIC CATALOGUE CUSTOMER'S AMENDMENTS ISSUE 07 — OCTOBER 2010

Page 1 of 8 pages

	Amendment
	NEW ITEMS
	Microcentrifuge, SCF2, Stuart Compact microcentrifuge for spinning 12 x 1.5/2.2ml tubes, or 4 x 0.2ml 8-tube strips using an accessory rotor.
	▲ Maximum speed 13500rpm (6000rpm - tube strip rotor)  ▲ 1 to 30 minute time  ▲ Safety lid interlock
	<ul> <li>▲ Pulse function</li> <li>▲ LCD readout of rpm or rcf and time</li> <li>▲ Quiet operation &lt;60dB</li> </ul>
	▲ Rapid rotor exchange without using tools
	Microcentrifuge, SCF2, Stuart As described. Supplied with 12x 1.5/2.2 microtube rotor (maximum rcf 12300xg). Overall 145 x 208 x 245mm H X W X D. Weight 4.4kg. For 220v 50/60Hz single phase supplies. CD170-15 SCF2 (Stuart SCF2)
	CD173-08 4 x 0.25ml 8-tube strip rotor (Stuart SCF2/2) (Maximum rcf 1850xg outer strip, 1542xg Inner strip.
	CD173-12 0.5ml tube adaptors for CD170-15 pk of 12 CD173-14 0.2ml tube adaptors for CD170-15 pk of 12
========	DELETIONS
=======================================	DELETIONS PP444-22 Cancelled, offer PP444-20
========	
========	PP444-22 Cancelled, offer PP444-20
	PP444-22 Cancelled, offer PP444-25  PP444-22 Cancelled, offer PP444-25  NEW ITEMS  Pasteur Pipettes Polyethylene. Disposable. With integral bulb. Graduated every 0.5ml. Length 155mm. Capactiy as indicated.  Sterile bulk wrapped in sleeves of 10. PP444-20 1ml, pack of 500
	PP444-22 Cancelled, offer PP444-20 PP444-22 Cancelled, offer PP444-25  NEW ITEMS  Pasteur Pipettes Polyethylene. Disposable. With integral bulb. Graduated every 0.5ml. Length 155mm. Capactiy as indicated. Sterile bulk wrapped in sleeves of 10.

# ISSUE 07 — OCTOBER 2010 Page 2 of 8 pages

### **AMENDMENTS**

TJ365-10 has a range of -20°C to +50 °C, not from -40 to +70 °C as shown in the catalogue.

TJ380-40 has a dimension of 68 x 37 x 18mm, rather than 55 x 25 x 105mm as shown in the catalogue.

All other details are unchanged.

### +-----

### **NEW ITEMS**

### **Hazard Warning Tape GHS-CLP**

In compliance with the U.N Globally Harmonised System (GHS) for classifying and labelling chemicals, under EU regulation EC1272/2008, Classification, Labelling and Packaging of hazardous substances and mixtures (CLP). Hard-wearing, self-adhesive vinyl in white printed with the pictogram as indicated and blank area to the side to allow the risk phase and hazard/contents information to be written or applied by the user. Pattern repeat 75mm with 50mm long blank area between symbols. Supplied singly in a containing roll of 66 metres long x 25mm wide on a 75mm core.

	CLP pictogram ref. GHS -	Hazard Class
SC250-12	GHS-01	Explosives, self-reactive substances, types A,B mixtures and organic peroxides
SC250-16	GHS-02	Flammable materials and flammable gas emitters, self-reactive or self-heating substances and mixtures, pyrophoric materials, organic peroxides
SC250-20	GHS-03	Oxidising substances
SC250-24	GHS-04	Compressed, liquefied, refrigerated, or dissolved gases
SC250-28	GHS-05	Corrosive agents
SC250-32	GHS-06	Acute toxicity (categories 1 to 3)
SC250-36	GHS-07	Acute toxicity (category 4), skin, eye or respiratory tract Irritation, sensitisation, organ toxicity, narcotic effects
SC230-40	GHS-08	Respiratory sensitizer, mutagenic/carcinogenic agent, reproductive/specific organ toxicity, aspiration harzard
SC250-44	GHS-09	Environmental toxicity

### **DELETION**

BS150-15 Cancelled, no replacement

## ISSUE 07 — OCTOBER 2010 Page 3 of 8 pages

### **DELETIONS**

WB200-40 Cancelled, replaced by WB200-42
WB205-45 Cancelled, replaced by WB200-42
WB210-38 Cancelled, replaced by WB210-44

WB215-43 Cancelled, replaced by WB210-44

#### **NEW ITEMS**

### Glassware Washers, Cub

Generally as previous WB200/WB205 models but with a backlit, graphics LCD displaying clock, calendar, timer and temperature readouts, alarm icons and troubleshooting menu, 10 programme controller and RS232 serial ports for printer and PC connection. With one each upper and lower level basket trays and standard spray arms. For 230V 50Hz single phase, 30A supplies.

WB200-42 Glassware Washer, Cub

Glassware Washer/Dryers, Cub

Generally as WB200/WB205 models, but with additional hot air drying unit.

WB210-44 Glassware Washer/Dryer, Cub Extra

Accessory Jet Racks

•	For	Hollow	For
	machines	spigots	levels
WB220-30	WB200/205	40 each	both
WB220-35	WB210/215	40 each	both

Baskets, Stainless steel mesh, wedge shape

For test tubes. For use on lower basket tray.

WB224-20 105mm high WB224-26 165mm high WB224-32 200mm high

Petri Dish Rack

Plastic coated stainless steel. For use on lower level only. For 31 dishes, 100mm diameter.

WB230-20 Petri dish rack

**Beaker/Bottle Support** 

Supports up to 28 items. Two supports can be used on each of the upper and lower basket trays.

WB232-10 Beaker/Bottle support

### **AMENDMENT**

FK500-37 is actually FK502-37

# Page 4 of 8 pages

· · · · · · · · · · · · · · · · · · ·	
	NEW ITEMS
	NEW TIENS
	Mixers
	MS224-09 Spare 2L stainless steel container and blending assembly with plastic lid
	MS224-13 Accessory stainless steel lid for MS224-09
	NIGHT-13 Accessory statiness seed in 101 Might-05
	AMENDMENT
	SB720-55 is now ref. LV54
	5B720-55 IS NOW Tell 12754
	All other details are unchanged.
	 +====================================
	AMENIDMENTEC
	AMENDMENTS
	WB890-series and WB950-series washbottles now have multilingual labelling in English, French,
	German and Spanish. They also include GHS hazard pictograms, NFPA diamond, risk phrases and
	safety information.
	All other details are unchanged.
	WB890-10 is now Azlon WGW602PML
	WB890-12 is now Azlon WGW603PML
	WB890-20 is now Azlon WGW600PML
	WB890-22 is now Azlon WGW601PML WB890-30 is new Azlon WGW530PML
	WB890-32 is new Azlon WGW536PML
	WB950-10 is now Azlon WGW531VTML WB950-12 is now Azlon WGW537VTML
	WB950-12 IS NOW AZION WGW53/VTML WB950-20 is now Azion WGW533VTML
	WB950-22 is now Azlon WGW539VTML
	WB950-30 is now Azlon WGW532VTML
	WB950-32 is now Azlon WGW538VTML WB950-40 is now Azlon WGW535VTML
	WB950-40 is now Azion WGW535 V IVIL WB950-42 is now Azion WGW541VTML
	WB950-50 is now Azlon WGW534VTML
	WB950-52 is now Azlon WGW540VTML
	DEV ETVON
	DELETION
	FZ545-35 Cancelled, offer replacement FZ545-15
	AMENDMENT
	FZ545-15 now also accommodates H400 hood
	All other late to an arrangement
	All other details are unchanged  ISSUE 07— OCOTBER 2010
	155UE V/— UCUIDER 2010

# Page 5 of 8 pages

	DELETION			
	BR102-50 Cancelled, no replacement			
	DELETION			
	DD4(0.40.C. H. I. 66. DD4(7.41			
	BD260-40 Cancelled, offer BD265-31			
=======		===========	========	==========
	DELETION			
	DC125 10 C			
	DS125-10 Cancelled, no replacement			
	DELETIONS			
	EV760 22 Concelled offer EV765 22			
	FK760-22 Cancelled, offer FK765-22 FK760-25 Cancelled, offer FK765-30			
	FK760-30 Cancelled, offer FK765-30			
	FK760-35 Cancelled, offer FK765-35			
	FK760-40 Cancelled, offer FK765-40			
	FK763-22 Cancelled, offer FK767-22			
	FK763-25 Cancelled, offer FK767-30			
	FK763-30 Cancelled, offer FK767-30			
	FK763-35 Cancelled, offer FK767-35			
	,			
	NEW ITEMS			
	NEW HENIS			
	Conical Flasks, Azlon, Screw Capped	l		
	Wide Neck. With PP screw cap and a	unnrovimato graduations. U	/ill also accent NS s	tonnore oe
	indicated. Autoclavable at 121°C. Su		m also accept 145 s	toppers as
		Capacity	Cap Size	Stopper accepted
	Polypropylene, transparent	ml	mm	NS
	1 orypropyrene, transparent			
	FK765-22	50	40	34/35
	FK765-25	100	40	34/35
	FK765-30	250	52	45/40
	FK765-35	500	52	45/40
	FK765-40	1000	52	45/40
	PMP (TPX) transparent		4.0	24/2-
	FK767-22	50	40	34/35
	FK767-25	100	40	34/35
	FK767-30	250 500	52 53	45/40
	FK767-35	500	52 52	45/40 45/40
	FK767-40	1000	52	45/40
	Accessory stoppers - see SW432/SW4	146 series		
		OCOTRED 2010		

### **CANCELLATION**

SG500-20 Cancelled, offer SG500-19

### **NEW ITEM**

Basic soil test kit

Generally as previous SG500-20 but with sufficient reagents to conduct soil tests. All other details are unchanged.

SG500-19 Basic soil test kit

#### **NEW ITEMS**

Filter Holders, Vacuum, Glass

Borosilicate glass with metal clamp, stainless steel membrane support screen, PTFE gasket and rubber stopper to fit accessory manifolds. Supplied singly.

	For membrane	Reservoir
	dia. mm	ml
FC335-05	25	15
FC335-09	47	300
FC335-13	90	1000

Filter Holders, Vacuum, Stainless Steel

All stainless steel construction for durability. Can be flame sterilised. With reservoir volume indicator mark and stopper to fit accessory manifolds. Supplied singly.

	For membrane	Reservoir
	dia. mm	ml
FC336-05	47	100
FC336-09	47	300
FC336-13	47	500

Accessory stainless steel lids for FC336 holders only. With vent. For funnels with reservoirs as indicated.

FC336-25	100ml	(167103-11)
FC336-29	300ml	
FC336-33	500ml	

### **Filter Manifolds**

Stainless steel. For use with FC335 and FC336 vacuum filter holders. With stopcock on each position. Supplied singly.

FC338-10	3 position
FC338-20	6 position

# T H E E N D

### PRICES FOR NEW ITEMS Page 7 of 8 pages

### **OCTOBER 2010**

CAT. NO.	PRICE £	Per	CODE
CD170-15	799.00	each	OS
CD173-08	100.00	pack	OS
CD173-12	70.00	pack	OS
CD173-14	70.00	pack	OS
FC335-05	113.00	each	00
FC335-09	143.00	each	00
FC335-13	430.00	each	00
FC336-05	332.00	each	00
FC336-09	355.00	each	00
FC336-13	408.00	each	00
FC336-25	51.00	each	00
FC336-29	51.00	each	00
FC336-33	56.00	each	00
FC338-10	1430.00	each	00
FC338-20	2740.00	each	00
FK765-22	8.50	each	AZ
FK765-25	9.52	each	AZ
FK765-30	9.92	each	AZ
FK765-35	12.83	each	AZ
FK765-40	16.87	each	AZ
FK767-22	10.64	each	AZ
FK767-25	10.93	each	AZ
FK767-30	12.98	each	AZ
FK767-35	15.94	each	AZ
FK767-40	19.00	each	AZ
MS224-09	245.00	each	ВО
MS224-13	35.00	each	ВО
PP444-20	16.30	pack	AL
PP444-25	16.30	pack	AL
SC250-12	4.10	roll	AO
SC250-16	4.10	roll	AO
SC250-20	4.10	roll	AO
SC250-24	4.10	roll	AO
SC250-28	4.10	roll	AO
SC250-32	4.10	roll	AO
SC250-36	4.10	roll	AO
SC250-40	4.10	roll	AO
SC250-44	4.10	roll	AO
SG500-19	18.00	each	FO

PRICES FOR NEW ITEMS

# Page 8 of 8 pages

WB200-42			5870.	00	each		ВО
WB210-44			8030.	00	each		ВО
	$\boldsymbol{T}$	H	$\boldsymbol{E}$	$\boldsymbol{E}$	N	D	