A J COPE & SON LTD The Oval London E2 9DU Telephone: +44 (0) 20 7729 2405 Facsimile: +44 (0) 20 7729 2657 Internet: www.ajcope.co.uk

ELECTRONIC CATALOGUE CUSTOMER'S AMENDMENTS ISSUE 7 — July 2012 Page 1 of 12 pages

Amendment
NEW ITEM <u>Bag Ties</u> Generally as BA095-10 but supplied in packs of 1000. BA095-15 Bag ties, pk.1000
CANCELLATIONS BD225-10 Cancelled, special order only BD225-20 Cancelled, special order only BD230-20 Cancelled, special order only NEW ITEMS Analytical Electronic Balances, HR-A/HR-AZ Series A&D Choice of autocalibration (HR-AZ) or external calibration using accessory weights (HR-A) models Weighing unit modes: g, mg, oz, ozt, ct, mo, dwt, gn, tl, tla, with %, counting, density and custom unit modes GLP-, GMP-, GCP- and ISO-compliant output Plastic draught shield with rotary sliding doors complies with FDA/HACCP "glass-free requirements" and unclips for easy cleaning and balance use in confined spaces Selectable weighing filters High contrast, reverse-backlit LCD readout

ISSUE 7 — July 2012 Page 2 of 12 pages

Mode	el	ext. cal. int. cal.	HR-100A HR-100AZ	HR-150A HR-150AZ	HR-250A HR-250AZ
Сара	city	g	102	152	252
Read	ability	g	0.0001	0.0001	0.0001
Repr	oducibility	±g	0.0001	0.0001	0.0001 (0 to 200g), 0.0002 (rest of range)
Linea	nrity	±g	0.0002	0.0002	0.0003
Pan s	ize	mm	90 dia.	90 dia.	90 dia.
Oper	ating temperature	°C	+5 to +40 at 8	5% RH or less —	- all models
Overa		W, mm	198	198	198
		D, mm	294	294	294
		H, mm	315	315	315
Weig	ht.	HR-A, kg	3.5	3.5	3.5
	,	HR-AZ, kg	3.9	3.9	3.9
BD2 BD2 BD2 BD2 BD2 BD2 BD2 BD2 BD2 BD2	rnal calibration mod25-40HR-100A25-45HR-150A25-50HR-250Aocalibration models (25-60HR-100AZ25-65HR-150AZ25-70HR-250AZ30-25USB adapte30-35LAN ethern00-55Calibration00-58Calibration00-62Calibration	(<u>internal weig</u> l r et adapter weight, F1, 50 weight, F1, 10	<u>nt)</u> g 0g		
	NSTATEMENT				
Bott	<u>les, Glass, Wincheste</u>	<u>er pattern</u>			
BR1		chester patter city, R6/31 nec		narrow mouth	n, without cap,
REI	NSTATEMENT				
Bott	les, Glass, Uncapped	<u>l</u>			
BR1	50-17 Bottles, med R3/22 neck,	-	g, amber glass	, 50ml capacit	у,

CANCELLATION BR154-35 Cancelled, offer BR154-40
REINSTATEMENT <u>Bottles, Powder</u> BR160-33 Bottles, powder, clear glass, wide mouth, without cap, 4oz capacity, R3/38 neck, pk.80
CANCELLATION BR404-35 Cancelled, offer BR404-33 REINSTATEMENTS <u>Bottles, Powder</u> BR400-33 Bottles, powder, clear glass, wide mouth with black plastic cap and PVDC liner, 4oz capacity, R3/38 neck, pk.80 BR404-33 Bottles, powder, clear glass, wide mouth with white polypropylene wadless cap, 4oz capacity, R3/38 neck, pk.80
NEW ITEM <u>Cleaner, Lipsol</u> A totally rinsable, phosphate- and enzyme-free, biodegradable surfactant with excellent cleaning performance on laboratory glassware, plasticware, rubber, ceramics and stainless steel. Mildly alkaline and therefore unsuitable for cleaning non-ferrous metals such as aluminium, zinc, or with polycarbonate plastics. Supplied in 5L plastic container. CL645-17 5 litres
 CANCELLATION FC372-34 Cancelled, offer FC378-35

CANCELLATIONS

LA470-50 Cancelled, offer LA477-45 LA470-55 Cancelled, no replacement LA470-58 Cancelled, no replacement LA470-62 Cancelled, no replacement LA473-03 Cancelled, no replacement LA473-07 Cancelled, no replacement LA473-11 Cancelled, no replacement LA473-15 Cancelled, no replacement LA473-23 Cancelled, no replacement LA473-27 Cancelled, no replacement

NEW ITEMS

Label printer, BMPTM51

Portable, thermal transfer label printer for labelling samples in accordance with GLP (Good laboratory Practice) protocols. Maximum label width printable 38.1mm. Print functions include text, bar-coding, serialisation, banner printing and multiple material configuration. Uses accessory drop-in cartridges containing a continuous reel of labels or tape up to 7.6m long. With dot matrix LCD, alphanumeric keyboard, label cartridge recognition system, USB connection/output port and type cutter. Overall 365 x 460 x 150mm. Weight 5.2kg. Supplied with hard carrying case, rechargeable Li-on battery, USB cable, control software, instructions, tape cartridge with black on white tape reel of 38.1mm width and a.c. adapter/charger requiring a 230V 50Hz single phase supplies.

LA477-45 BMPTM51

Accessory tapes

Supplied in drop-in cartridges containing reels of labels as indicated. All print black text on white background.

	Label size, mm	Material range	Temp. °C	For labelling	Labels per cartridge
LA478-04	9.53 dia.	Polyester	-196 to +121*	0.5 to 1.5ml microtubes	240
LA478-08	25.4 x 9.53/ 9.53 dia.	Nylon	-196 to +121	1.0 to 2.0ml microtubes	240
LA478-12	12.7 dia.	Nylon	-196 to +121	1.5 to 2.0ml microtubes	240
LA478-16	38.1 x 98.42	Polyester	-196 to +121	50ml tubes	80
LA478-20	25.4 x 48.26	Polyester	-196 to +121	12 to 24ml tubes	100
LA478-24	25.4 x 48.26	Polyester	-196 to +121*	50ml tubes	100
LA478-28	38.1 x 31.75	Polyester	-196 to +121*	12 to 24ml tubes	160
LA478-32	6.35 x 48.26	Polyester	-80 to +121	Microplates	160
LA478-36	19.08 x 12.70	Vinyl	-80 to ambient	Tissue cassettes	220
LA478-40	22.86 x 22.86	Polyester	-80 to +121	Microscope slides	280

*Not recommended for use in boiling baths over medium to long periods.

CANCELLATIONS

MW605-40 Cancelled, offer MW605-50 MW605-45 Cancelled, offer MW605-55 MW606-15 Cancelled, offer MW606-25 MW606-20 Cancelled, offer MW606-30 MW606-40 Cancelled, offer MW606-50 MW606-45 Cancelled, offer MW606-55 MW607-26 Cancelled, no replacement MW607-29 Cancelled, no replacement MW607-50 Cancelled, no replacement

NEW ITEMS

<u>Personal Thermal Cyclers ³Prime Series</u> <u>Techne</u>

Self-contained, personal thermal cyclers with numerous applications in molecular biology and diagnostics laboratories. A choice of models is offered with fitted blocks either for 0.2ml or 0.5ml microtubes in strips or plates. A full colour touchscreen provides control and display of programme status, with templates allowing easy programming of even complex procedures. Up to 1000 x 3-segment programmes can be stored with incremental/decremental timing and temperature elements. Programme naming and password protection are also available. A heated lid is provided with temperature adjustment (100 to 115°C) and heat on/off switch. The heated lid for ³PrimeX and ³PrimeG models also adjusts to accommodate a variety of consumables. Additionally, ³PrimeG models are convenient, small-scale, gradient thermocyclers which will also hold half a 96-well plate in an 8 column x 6 row configuration. Annealing temperatures can be optimised over a 1 to 14°C gradient between 30°C and 80°C. All ³Prime models include a USB output to transfer programmes and temperature logs between instruments or the user's computer. For 100 to 230V 50/60Hz single phase supplies.

Model		³ Prime ³ PrimeX	³ PrimeG
Temp. range, block	°C	+4 to ± 100 — all models	Timeo
gradient	°Č		+30 to +80
Block uniformity (at 55°C)	°Č	<±0.3 — all models	
Temperature accuracy (at 55°C)	°Č	$<\pm 0.25$ — all models	
Min./Max. gradient	°Č		1/14
Max. No. of programmes [†]	-	1000 built-in — all models	_,
(3-segment)			
Run end time calculation		Yes — all models	
Selectable heated lid temp. [#]	°C	100 to 115 — all models	
End of programme alarm		Yes — all models	
Manual pause and stop facility		Yes — all models	
Incremental/decremental temp.		Yes — all models	
control			
Programme ramp rate,		up to 3.0 — all models	
°C per second		*	
Peltier elements per block		4 — all models	
Data/network interface		USB — all models	
Dimensions D x W x H	mm	350 x 210 x 180 — all models	
Weight, net	kg	6 — all models	
	C		

[†]The maximum number of programmes available depends on the number of segments required per programme.

[#]Heated lid will only come on if block temperature is set above 35°C.

Personal Thermal Cyclers, 3Prime, Techne

MW605-50 ³Prime with 24 x 0.2ml microtube block MW605-55 ³Prime with 18 x 0.5ml microtube block

MW606-25 ³PrimeX with 48 x 0.2ml microtube block MW606-30 ³PrimeX with 30 x 0.5ml microtube block

MW606-50 ³PrimeG with 48 x 0.2ml microtube block MW606-55 ³PrimeG with 30 x 0.5ml microtube block

Accessory

³Prime Gradient Upgrade

Allows ³PrimeX models only to be upgraded to support gradient operation.

MW606-95 ³Prime gradient upgrade

Note: alternative replacement interchangeable blocks can be fitted by returning the instrument to ourselves for factory service.

CANCELLATIONS

MW630-27 Cancelled, offer MW630-45 MW630-33 Cancelled, offer MW630-40 MW633-05 Cancelled, no replacement MW633-08 Cancelled, no replacement MW633-16 Cancelled, no replacement MW633-20 Cancelled, no replacement MW635-26 Cancelled, offer MW635-45 MW635-44 Cancelled, offer MW635-35 MW637-05 Cancelled, no replacement MW637-08 Cancelled, no replacement MW637-16 Cancelled, no replacement MW637-20 Cancelled, no replacement

NEW ITEMS

<u>Programmable Thermal Cyclers Prime Series</u> <u>Techne</u>

High performance, thermal cyclers designed for laboratories requiring a large sample throughput. A fully interchangeable block system allows exchange without needing tools and offers enhanced versatility. A large, 5.7" VGA colour touchscreen with intuitive software programming and graphical display ensures user-friendly control and operation. Up to 1000 x 3-segment programmes can be stored with incremental/decremental timing and temperature elements. 8 Peltier modules per block provide rapid heating and cooling over the range $+4^{\circ}$ C to $+100^{\circ}$ C, with block uniformity $<\pm0.3^{\circ}$ C at 55°C. A heated lid is provided with adjustable sealing pressure to accommodate a variety of consumables and a USB port allows networking with other cyclers or data transfer to the user's computer. Additionally, the PrimeG model allows gradient cycling with annealing temperatures optimisable over a 1 to 29°C range (1 to 19°C when using the 384-well block) and between 30°C and 80°C. For 100 to 230V 50/60Hz single phase supplies.

Model	0.0	Prime	PrimeG
Temp. range, block	°C	+4 to +100 — both mode	els
gradient	°C	_	+30 to +80
Block uniformity (at 55°C)	°C	<±0.3 — both models	
Temperature accuracy (at 55°C)	°C	<±0.25 — both models	
Min./Max. gradient	°C	_	1/29 (19 - 384 well
-			plate block)
Max. No. of programmes [†]		1000 built-in — both ma	dels
(3-segment)			
Selectable heated lid temp. [#]	°C	100 to 115 — both mode	ls
Programme ramp rate,		up to 3.4 — both models	
°C per second			
Peltier elements per block		8 — both models	
Data/network interface		USB — both models	
Dimensions D x W x H	mm	420 x 240 x 240 — both	models
Weight, net	kg	9.4 — both models	

[†]The maximum number of programmes available depends on the number of segments required per programme. [#]Heated lid will only come on if block temperature is set above 35°C.

<u>Programma</u>	ble Thermal Cycler	rs, Techne	
Prime, non-	<u>gradient</u>		
	Prime for 96 x 0.2r		
	Prime for 60 x 0.5r Prime for 384-well		
PrimeG, gra			
rimeG, gra	<u>idient</u>		
	PrimeG for 96 x 0. PrimeG for 60 x 0.		
	PrimeG for 384-we		
Spares and a	accessories		
MW637-40	Accessory block fo	r 96 x 0.2ml microtubes	
	v	or 60 x 0.5ml microtubes	
	•	or 384-well microplates	
MW637-60	Prime gradient upg gradient models	grade. Converts Prime cyclers to Prime	
	grautent models		
 CANCELLA	 ATION		
PD100-65 C	ancelled, no replace	ement	
 			===
NEW ITEM	IS		
Splash head	adapters, rotary ev	aporator	
For use with	n rotary evaporators	s to prevent "bumping" or foaming solvents from lable with, or without, sintered glass (porosity 1) fi	lter.
		one size	
	Flask	Condenser	
<u>Without filt</u>	<u>er</u>		
QSH15/19	19/26	29/32	
QSH15/24	24/29	29/32	
QSH15/29	29/32	29/32	
<u>With filter</u>			
QSH16/19	19/26	29/32	
QSH16/24	24/29	29/32	
QSH16/29	29/32	29/32	

	NEW ITE	MS				
	Pharmacy	Refrigerators				
	<u>Liebherr</u>					
		e and medicin as. Internal te		harmaceutical dispe 5°C ±3K.	nsaries, surger	ies and other
	 Digital Safety Volt-fr Power Access Seamle Lockat 	display with r thermostat to ee and RS485 failure alarm port for indep ss plastic inne ole solid or gla	nin/max men ensure tempe serial interfa indicator pendent temp er liner for ea ass door mode	erature does not fall l ce erature probe sy cleaning	below 2°C	
			0	, door lock and adju Hz single phase supp		esh wire
		Model	Gross cap. litres	H x W x D, mm	Shelves supplied	Internal light
	Solid door m	odels				
	RG285-10 RG285-30	MkUv1610 MkUv3910	141 360	820 x 600 x 615 1840 x 600 x 615	3 5	No No
	<u>Glass door m</u>	odels				
	RG285-60 RG285-80	MkUv1612 MkUv3912	141 360	820 x 600 x 615 1840 x 600 x 615	3 5	Yes Yes
:=======	CANCELI	ATIONS				
		Cancelled, off Cancelled, off				
	NEW ITE	MS				
	<u>Steristopp</u> Generally		W536-series b	ut with reduced pacl	x sizes.	
		To fit tube i	.d. mm Pac	k qty.		
	SW536-21 SW536-39	10 15	100 500			

==

	CANCELLATI	ON	
	SX240-35 Canc	elled, offer SX2	40-45
	NEW ITEM		
	case and lid are equipped with s 79mm high.	pylene case wit gridded to hold	h clear polycarbonate, hinged cover. The floor of the l 15, removable, full depth inserts in place. The case is atches and hinged, carrying handle. Overall 413 x 330 x
	5A240-45 Ca	TTy-Lite Case	
============	NEW ITEMS		
	<u>Culture Flasks,</u>		
	Generally as TI shaking of cultu		cal flasks but with baffles to increase agitation during
		Ref.	Capacity, ml
	TP585-05	1134/08	250
	TP585-10	1134/10	500
	TP585-15	1134/12	1000
	TP585-20	1134/14	2000
	Fernbach Flask		
	Generally as TI	2585-series but	Fernbach shape with baffles.
		Ref.	Capacity, ml
	TP585-80	1138/10	2800
		TH	E END

PRICES FOR NEW ITEMS Page 11 of 12 Pages

uly 2012	Page 11 of 12 Pag	ges	
CAT. NO.	PRICE £	Per	CODE
BA095-15	18.95	pack	
BD225-40	765.00	each	RW
BD225-40 BD225-45			
	825.00	each	RW
BD225-50	980.00	each	RW
BD225-60	930.00	each	RW
BD225-65	975.00	each	RW
BD225-70	1140.00	each	RW
BD230-25	83.00	each	RW
BD230-35	438.00	each	RW
BR120-50	24.60	pack	RO
BR150-17	20.90	pack	RO
BR160-33	22.00	pack	RO
BR400-33	33.15	pack	RO
BR404-33	29.75	pack	RO
CL645-17	30.90	each	AZ
	25.74	each per 2	
LA477-45	459.00	each	FO
LA478-04	62.30	pack	FO
LA478-08	32.90	pack	FO
LA478-12	32.15	pack	FO
LA478-16	45.40	pack	FO
LA478-20	70.30	pack	FO
LA478-24	26.15	pack	FO
LA478-28	41.25	pack	FO
LA478-32	28.95	pack	FO
LA478-36	26.55	pack	FO
LA478-40	49.20	pack	FO
MW605-50	2195.00	each	ED
MW605-55	2195.00	each	ED
MW606-25	2795.00	each	ED
MW606-30	2795.00	each	ED
MW606-50	3195.00	each	ED
MW606-55	3195.00	each	ED
MW606-95	400.00	each	ED ED
MW630-40	3795.00	each	ED ED
MW630-40 MW630-45	3795.00	each	ED ED
MW630-50	3795.00	each	ED
MW635-35	4395.00	each	ED
MW635-45	4395.00	each	ED
MW635-50	4395.00	each	ED

PRICES FOR NEW ITEMS Page 12 of 12 Pages

MW637-40	2125.00	each	E
MW637-45	2250.00	each	E
MW637-50	2250.00	each	E
MW637-60	600.00	each	E
QSH15/19	61.80	each	C
QSH15/24	63.86	each	C
QSH15/29	65.92	each	C
QSH16/19	67.98	each	C
QSH16/24	70.04	each	C
QSH16/29	70.00	each	C
RG285-10	1195.00	each	Μ
RG285-30	1565.00	each	Μ
RG285-60	1285.00	each	M
RG285-80	1725.00	each	M
SW536-21	39.73	pack	A
SW536-39	23.00	pack	A
SX240-45	66.60	each	R
TP585-05	21.22	each	C
	10.61	each per 5	
TP585-10	25.46	each	C
	12.73	each per 5	
TP585-15	47.75	each	C
	31.83	each per 2	
TP585-20	42.44	each	C
TP585-80	200.85	each	C
	133.90	each per 2	
T	HE E	CND	